
SUBLIMADOR MONODOSIS SBL003

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OXALIC ACID VAPORIZER SBL003

DESCRIPTION

- 1 Temperature regulator
- 2 Dispenser
 - 2a Push
- 3 Rear support rod
- 4 Dispenser receiving cap
- 5 Oxalic acid vaporizer vesse
- 6 Ejector spout
- 7 Thread sleeve for securing the rear support rod

GENERAL RECOMMENDATIONS

Sublimation is the transition of a substance from the solid to a gaseous state without passing through a liquid state. When heated, oxalic acid decomposes into toxic gases such as carbon monoxide and carbon dioxide.

- Not suitable for persons under the age of 18.
- Before carrying out any operation with this machine you must be fully aware of how it works.
- Before operating you must wear: disposable coverall suitable for acids, full face gas mask, rubber gloves suitable for acids, boots suitable for acids.

**Filters 3M™ 6059
Recommended
multigas filter**

**Electrical
HAZARD**

**Burn
HAZARD**

**WARNING
corrosive**

**Inhalation
HAZARD**

**Respirators
must be worn
in the area**

INSTRUCTIONS FOR USE

Connect the machine to the power supply, set the temperature (1) and wait for it to reach the desired temperature. Fill the dispenser (2) with the desired amount of oxalic acid. Insert the dispenser into the receiving cap (4) and gently press the push (2a). Insert the ejector spout (6) into the hive entrance and wait until the gas dispensed by the machine is depleted.

TEMPERATURE REGULATION (FIG.1)

Once the oxalic acid vaporizer is connected, the current protector turns on, checks the voltage (A) [10 sec.], the thermostat turns on (B) [3 sec.] and the heating element begins to heat up*. Working temperature: 260°C**.

To change the thermostat temperature:

1° Press SET(C) once.

2° Press -(D) or +(E) until the desired temperature is set.

3° Press SET(C) again.

4° The machine will adjust the heating element to the desired temperature.

- ¿What is the reason for modifying the temperature?:

The ambient working temperature. In environmental conditions below 6°C, it is recommended to raise the temperature by 10°C-15°C.

- ¿How do I know it is sublimating properly?:

Place a wooden board approximately 10 cm from the ejector spout, sublimate and observe that the product should stain the board white.

(*) The temperature indicated does NOT refer to the sublimation temperature, but rather to the resistance of the sublimation vessel.

(**) It maintains a constant temperature in the sublimation vessel, ensuring proper sublimation throughout the treatment.

REGULATION OF THE DISPENSER (FIG.2)

Although the dispenser (F) is initially set to load 4 grams, you can adjust the dosage using the adjustment nuts (K) [1-6 grams maximum].

- It is recommended not to wet the tip of the dispenser with water during sublimation. If a small crust forms, remove it with a clean, damp cloth.

- If oxalic acid gets between the sealing gasket (J) and the pusher (L), gently tap the dispenser to remove it and avoid reducing its load capacity.

- After prolonged use of the machine, check the condition of the recoil spring (I) to ensure that it continues to provide adequate resistance to the dispenser push button (G), allowing for optimal oxalic acid loading.

- You can modify the travel of the dispenser push button using the travel nuts (H).

REPLACEMENT OF THE DISPENSER RECEIVING CAP & CLEANING OF OXALIC ACID VAPORIZER VESSEL (FIG.3)*

Remove the retaining screws (M) [Allen no. 5] to replace the dispenser receiver cap (N) or access the oxalic acid vaporizer vessel (O). If the machine loses efficiency, check the condition of the oxalic acid vaporizer vessel. If there is a crust stuck to the inside. To remove it, maintain the working temperature in the machine, pour a little water mixed with universal degreaser (ratio 2:1) into the oxalic acid vaporizer vessel, quickly place the dispenser and allow the mixture to evaporate. Do this several times to remove all traces. Finish cleaning by repeating the process two or three times with water only to leave the oxalic acid vaporizer vessel clean.

(*) Do not tighten the screws too much; a moderate tightening is sufficient, as the temperature adjusts the parts.

CLEANING THE OXALIC ACID VAPORIZER

Clean the machine at the end of each use. Disconnect the unit from the power supply and, while the machine is still hot, pour water into the opening in the receiving cap and quickly replace the dispenser, thereby increasing the pressure inside the sublimation vessel. Allow the water to evaporate through the oxalic acid vaporizer's ejector spout. Repeat the process three or four times. This will resolve most blockages. **WARNING:** high pressure and hot steam.

IMPORTANT NOTES

- If the machine is connected to auxiliary power equipment (generator), this must be running before connecting the sublimator, in continuous operation mode, and must have overvoltage protection.
- The machine has been tested beforehand to make the necessary adjustments and ensure optimal performance.
- Wait for the machine to reach operating temperature. Once sublimation begins, the temperature drops, the internal regulator readjusts the temperature, and the machine continues to operate normally, without the need to stop sublimation to increase the temperature.
- DO NOT operate in damp or rainy conditions. Risk of electric shock.
- DO NOT hang the machine from the spout without using the rear support rods, as this may cause the ejector spout to break.
- DO NOT connect and leave the machine to heat up on a heat-sensitive surface. Risk of fire.
- DO NOT leave the machine connected without supervision. Risk of fire.

WARRANTY AND REGISTRATION

- Any tampering with the machine without the manufacturer's written authorisation shall be understood as termination of the warranty period.
- Misuse of the machine may void the warranty.
- Scan the following QR code to register your warranty:

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